

Lessons Learned Regarding Delivery of Aerosolized Medicines in Patients With COVID-19

Original article: Ari A, Scott JB. Lessons learned about aerosol drug delivery in the era of COVID-19. Respiratory Care NetWork. Available at: <https://www.chestnet.org/Topic-Collections/COVID-19/COVID-in-Focus/Lessons-Learned-About-Aerosol-Drug-Delivery-in-the-Era-of-COVID-19>. Published July 30, 2021. Accessed September 30, 2021.



Over the course of the COVID-19 pandemic, multiple lessons have been learned regarding the delivery of aerosolized medicines



Lesson 1

The risk of device contamination and viral transmission is not the same across all devices

Contamination may be less with a mesh nebulizer versus a jet nebulizer

Mesh nebulizer

- Medication reservoir is separate from the patient interface
- Electric/battery powered
- Small residual volume

Jet nebulizer

- Open reservoir that is placed below the circuit → presents a risk of contamination via infected condensate
- Gas powered → may increase exhaled aerosol dispersion
- Large residual volume → may be hospitable to pathogens



Lesson 2

An aerosol device should be selected based on the patient's clinical status

Use a **mesh nebulizer** or a valved T-piece with a jet nebulizer in patients receiving HFNC therapy, NIV, or MV, to **keep the ventilator circuit intact**



Lesson 3

Aerosol delivery via a HFNC should be considered

Use a **mesh nebulizer** placed prior to the humidifier



Lesson 4

In COVID-19 patients, device and interface selection are of equal importance

Use a mouthpiece for aerosol delivery in spontaneously breathing patients



Lesson 5

Adopt good infection control and prevention measures to reduce exhaled aerosol dispersion to the environment

- Attach a filter to the expiratory outlet of nebulizers and ventilators
- Place a surgical mask over a HFNC to decrease dispersion of aerosols

Lessons Learned Regarding Delivery of Aerosolized Medicines in Patients With COVID-19

Original article: Ari A, Scott JB. Lessons learned about aerosol drug delivery in the era of COVID-19. Respiratory Care NetWork. Available at: <https://www.chestnet.org/Topic-Collections/COVID-19/COVID-in-Focus/Lessons-Learned-About-Aerosol-Drug-Delivery-in-the-Era-of-COVID-19>. Published July 30, 2021. Accessed September 30, 2021.

Spontaneous breathing	HFNC	Mechanical ventilation	Tracheostomy
Use aerosol drug delivery only if necessary	Use a HFNC for aerosol delivery before severe hypoxemic respiratory failure develops	To keep the ventilator circuit intact and to prevent viral transmission, do not use a pMDI	To prevent viral transmission, do not use a pMDI placed directly to the tracheostomy tube
Use an inhaler in patients who can optimally use one	Use a mesh nebulizer placed before the humidifier and with the reservoir cap kept closed	Use a mesh nebulizer placed before the humidifier and with the reservoir cap kept closed	Use a mesh nebulizer or a jet nebulizer with a valved T-piece, attached to the high-flow oxygen delivery system
Use a nebulizer with a filter placed on the expiratory port in patients who cannot optimally use an inhaler or who require drugs not available in an inhaler	Ensure that the prongs of the HFNC fit well and are not loose	Use a jet nebulizer with a valved T-piece if a mesh nebulizer is not available	Do not deliver aerosol therapy via the tracheostomy mask
Use a mouthpiece (not a facemask)	Place a surgical mask over the HFNC prior to aerosol delivery	Add a HEPA filter to the ventilator's exhalation arm and ensure the ventilator circuit remains intact	Deliver aerosol using an unassisted instead of an assisted technique unless a filter has been placed on the expiratory port of the manual resuscitation bag

For all patients during aerosol drug delivery

Isolate patients, administer aerosols in a negative-pressure room, wear appropriate PPE, and observe aseptic technique when preparing, cleaning, and maintaining all devices

Want to know more? Scan or click the QR code

