

OxyMask™




OxyMulti-Mask™ Instructions



Proper positioning and fit of the OxyMulti-Mask™ is important.

An appropriately sized OxyMulti-Mask™ will fit below the eyes and cup the underside of the chin. **Center** the pin and diffuser between the nose and above the upper lip area.

- 1 **Remove the mask from the packaging, recycle or discard the plastic insert.**
- 2 Twist swivel adapter until it rotates freely.
- 3 **For oxygen therapy**, ensure connector is attached to the swivel elbow of the mask and attach the oxygen tubing to the tapered end. Attach the universal connector at the end of the mask tubing to an oxygen source.
- 4 Place the mask on the patient's face and position elastic strap at a comfortable tension below the ears or at the crown of the head. Adjust tubing at the swivel to relieve any twisting.
- 5 **Adjust the oxygen flow from 1 to 15+ lpm to achieve prescribed patient oxygen saturation or oxygen flow rate.**
- 6 Allow patient to stabilize on oxygen flow and assess. Adjust oxygen flow rate as necessary to achieve or maintain the prescribed oxygen saturation.
- 7 **For medicated aerosol therapy**, remove the connector from the mask's swivel elbow, attach a nebulizer set up, push up and twist to secure. Attach the oxygen tubing to the flow meter and adjust the flow between 6 to 8 litres per minute, or as per nebulizer manufacturer's specifications.



OxyMulti-Mask™ Flow Rates

Actual FI_{O_2} will vary. Breath rate, depth of breathing and distance of the diffuser from the face will affect FI_{O_2} .

Actual FI_{O_2} required should be prescribed in conjunction with the use of a pulse oximeter.

Flowrate	FI_{O_2}
1 LPM	23% - 26%
2 LPM	27% - 33%
3 LPM	32% - 48%
4 LPM	45% - 57%
5 LPM	56% - 68%
7 LPM	64% - 74%
10 LPM	70% - 77%
12 LPM	70% - 81%
≥15LPM	72% - 83%



- Respiratory
- Surgical
- Blades & Scalpels
- Dynamic Tissue Systems

Learn more at: thebetteroxygenmask.com
oxygenmask@southmedic.com